DEUTSCH

ENGLISH

FRANÇAIS

NEDERLANDS

ESPAÑOL

ITALIANO

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# VM 9100RM III


Branding: **ANDUMEDIC® 3**


Gebrauchsanweisung  
User Manual  
Manuel d'utilisation  
Gebruiksaanwijzing  
Manual de instrucciones  
Istruzioni per l'uso  
инструкция по эксплуатации

Made in Germany

CE 0633

Aufbauanleitung  
Mounting instructions  
Instructions de montage  
Montagehandleiding  
Instrucciones de montaje  
Istruzioni per il montaggio  
Руководство по сборке

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hhp  
home health products  
[www.hhp-international.com](http://www.hhp-international.com)

*Dear customer*

*Thank you for purchasing this therapy system and for your trust in our company.*

*To make sure that you can take full advantage of all the benefits of our therapy system, we have compiled all the information which is important for optimum use in these operating instructions.*

**The operating instructions also contain important information on SAFETY and TECHNOLOGY. Please read them thoroughly before commissioning and store them in a safe location.**

*For this little effort, you will be rewarded by a long service life and the error-free function of our therapy system.*

## **Contents**


Assembly instructions.....	1/2
Test marks and markings.....	4/5
Technical data.....	5
Intended use.....	6
Accessories.....	6
Commissioning.....	6
Equipment.....	7
Technology and design.....	7
Infrared heat.....	8
Manual control.....	9
Buttons and indicators.....	10/11
Body zones.....	11
Pre-programmed intensity levels in the default settings.....	12
Program P1 – P20.....	12/13
Safety instructions.....	14
User instructions.....	15
Environmental protection.....	16
Power down.....	16
Cleaning, care and disinfection.....	17
Troubleshooting.....	18
Guidelines on electromagnetic compatibility – EMC guidelines.....	19-23
Addresses.....	23

## Test marks and markings of the VM 9100RM III therapy system


---

# CE 0633

Medical device conformity marking  
Notified body Berlin Cert GmbH

Max. operating time (therapy time)

After maximum operating period of 60 min  
(Therapy time) observe a break of 10  
minutes.

Applied part type BF  
Bodyfloating (unearthed)

Year of manufacture

**DC** Direct current

Electronic waste – Do not dispose  
of in household waste

Observe the operating instructions!  
(white symbol on blue  
background)

Manufacturer

Keep away from moisture and  
humidity “umbrella”

**IP21**

Protected against solid  
foreign objects  $> \varnothing 12.5 \text{ mm}$ ,  
protected against dripping water

Caution

## Test marks and markings: power supply unit (see technical data, page 5)

---

Direct voltage

Polarity, direct voltage power supply unit.  
The polarity must match the device used.

TÜV Rheinland standards testing  
acc. to EN60601-1/IEC 60601-1

Electronic waste – Do not  
dispose of in household waste

US test mark (Recognized Component  
Mark), tested for Canadian and US  
markets

Communauté européenne  
(European community).  
Directives for electrical and  
electronic devices.

## Test marks and markings: power supply unit (see technical data, page 5)


---

Do not use outdoors

FCC is the US certification mark for electromagnetic compatibility by the Federal Communications Commission.

Double insulation, contact protection, protection class II

## Technical data


---

Dimensions: spread out: approx. 1960 x 750 x 85 (max. 350) mm (L x W x H)  
folded: approx. 1130 x 750 x 175 mm (L x W x H)

Weight: approx. 19 kg

Rated voltage: 12 V DC

Rated power: approx. 80 W

Power supply unit: EM 11202E  
prim. 100-240 V AC, 2.0-1.0 A, 50-60 Hz  
sec. 12 V  8.33 A, max. 100 W  
Manufacturer: EDAC / METEK GmbH

IP rating: IP21 (for therapy system incl. accessories)

Maximum admissible body weight of the user: 120 kg

Maximum achievable surface temperature of the device during application: 44 °C

Firmware version: V1.1 P513P1\_3

Service life: 10 years with an average operating time of max. 2 hours per day.

Test marks and markings:

**CE 0633**

Manufacturer/  
service address: METEK GmbH  
Stammweg 8  
37327 Leinefelde-Worbis  
Germany

## Intended use

---

The therapy system is intended for the treatment of muscles and joints, for the vitalization of tissue which accumulates liquid and for regeneration.

The patient is an intended operator.

## Accessories:

---

- Power supply unit, (see technical data, page 5), fixedly attached to the device
- Wedge-shaped cushion (health cushion)
- Infrared – neck cushion (IR-NK)
- Infrared – back cushion (IR-RÜ)
- Infrared – calf cushion (IR-WA)
- Andullation abdominal belt (AG)
- Power – andullator (POW)
- Support wedge (Heart – weighing cushion)

## Commissioning

---

- First, check if all packaging is free of defects/damage. Carefully remove the therapy system and the accessories from the boxes or polybags and check whether all components are intact.
- Unfold the system and place it on an even surface.
- Connect plugs and accessories according to the assembly instructions (see page 1 and 2).
- Plug the mains supply line of the power supply unit in the socket.
- The flashing LED indicator next to the OFF button  on the manual control indicates that the device is ready for operation.
- Lie down on the therapy system so that the logo is located at the foot.
- Now pre-select a duration of 15 or 30 minutes using the watch button .
- After that, press the program button P/ON.  The program mode is activated.

Now, use the + / - buttons  to select the program.

The selected program number is shown on the display.   
After approx. 1.5 seconds, the therapy system begins to massage.

- The INFRARED HEAT, PROGRAMS, INTENSITY and TIME functions are now freely selectable for you to configure your individual massage program.
- This device is not intended for continuous operation. After the maximum therapy time (60 minutes), there must be a minimum break of 10 minutes.

### Note:

To switch off the therapy system completely, pull out the mains plug with all poles. This plug must be easily accessible and unobstructed. Do not pull the plug by pulling at the cable. Only use a properly installed socket.

If you need help commissioning or using the system or would like to report on unexpected situations / events during operation, please contact the manufacturer/service address indicated at the end of these operating instructions.

## **Equipment**

---

- State-of-the-art comfortable manual control with background lighting and sleep mode
- Neck, back and calf cushion with infrared.
- Andullation belt, heart balance cushion and integrated arm massage
- Unicoloured or bicoloured artificial leather, tested for irritations and sensitisation as well as cytotoxicity with the following functional elements:
  - Support cushion for head and neck
  - Reinforcing ribs to intensify the massage
  - Nub surface for additional accentuation of the massage in the foot, leg, back, belly and neck area.
  - Adjustable infrared pillow
  - Loop to hold the manual control
  - Connectible infrared heat

## **Technology and design**

---

- The therapy system stands out due to sophisticated characteristics such as state-of-the-art technology, elegant design and high functionality. Many supposedly minor details add up to its superior quality.
- Regarding the design of the product, safety was always the main concern. This is why the therapy system is equipped with the following safety characteristics:
  - › The motors of the vibration units and the infrared heat/heating system are operated by safety extra-low voltage
  - › Fuse for protection in case of defective or sluggish motors
  - › Foam, CFC-free
  - › Artificial leather in heat-resistant design
  - › Sleep mode
  - › Timer operation
  - › connected power supply unit
- The core of the therapy system is a microprocessor control which assigns your input commands via micro-buttons to the corresponding function.
- The vibration units are embedded in foam which is enclosed by a cover based on the different zones of the body.
- The modular design makes an important contribution to environmental protection. Many materials can be retrieved and recycled in a recycling process.
- For the following components of the VM 9100RM III therapy system, cables are available/connected:
  - › Manual control
  - › Power andullator
  - › Infrared neck cushion
  - › Infrared back cushion
  - › Infrared calf cushion
  - › Abdominal andullation belt
  - › Power supply

## Infrared heat

---

**Caution: Select the heat level according to sensitivity and skin type**


- The radiation penetrates the skin deeply. It is heated so that the blood vessels are dilated and the circulation is stimulated.
- The heat effect facilitates the opening of pores.
- The infrared heat helps protect the health of your skin since it has a stimulating effect on the blood circulation.

**IMPORTANT!**

- Do not expose your eyes or other sensitive body parts to infrared heat. Improper use can lead to burns.
- Persons who receive medical treatment or suffer from inflammations, swelling or rash should consult a physician before using infrared heat.


## Manual control

## Buttons and indicators

### *OFF button:*

By operating the “OFF” button, you stop all functions. The readiness for operation is still indicated by the flashing LED next to the “OFF” button.

### *+ / - button*

Button for individual settings + / -

### *Watch button*

By pushing this button, you select the pre-programmed massage duration of 15 or 30 min **before selecting a program**. During the massage, you activate the TIME MODE and can individually set the time (5-60 min) using the +/- buttons.

### *Display*

Program, time and intensity display. After a change of the setting, the display always returns to the current program after 5-10 sec. When the set program is completed, -End- is displayed.

### *PROGRAM button*

By pushing this button, the PROGRAM MODE is activated. (See also program list P1 – P20) You can now set the desired program using the +/- buttons. It is, of course, possible to change between **P1** and **P20** while a program is running using the +/- button. The selected program is shown on the display.

**NOTE:** The program sequence is designed to let approx. 1.5 seconds pass when a program P1 – P20 has been selected before the massage starts.

### *TIME*

LED display – Time

### *INT*

LED display – Intensity

### *PROG*

LED display – Program

### *Zones button 1-5*

By pushing these buttons, the INTENSITY MODE of the individual zones 1-5 is activated (only possible after program selection or during the massage). Now, you can select the desired intensity using the +/- buttons. 7 intensity levels are available. The set intensity level is shown on the display for approx. 5-10 s. The LEDs 1-5 are assigned to massage zones 1-5. The LED 1-5 of the message zone which is currently active is lit. The +/- buttons can also be used to deactivate individual massage zones. The corresponding LED is then switched off and “OFF” is shown on the display.

### *Zones complete button*

By pushing this button, the INTENSITY MODE of zones 2-5 is activated. Using the +/- buttons, the desired massage zones 2-5 can be set. This setting can be made prior to program selection and during the massage. 7 intensity levels are available. The selected intensity level is shown as “mit Ic 1-7” on the display for approx. 5-10 seconds.

H 1 button Using this button, the infrared heat of the back and neck cushion can be set to four levels immediately after the system is ready for operation and in any program:

- Push 1 time → level 1– low heating power → H1 LEDs flash slowly
- Push 2 times →level 2 – medium heating power → H1 LEDs flash quickly
- Push 3 times →level 3 – maximum heating power → H1 LEDs are permanently lit
- Push 4 times → infrared program is activated → H1 LED double flashing interval

H2 button Using this button, the infrared heat of the calf cushion can be set to four levels immediately after the system is ready for operation and in any program:

- Push 1 time → level 1– low heating power → H1 LED flashes slowly
- Push 2 times →level 2 – medium heating power → H1 LED flashes quickly
- Push 3 times →level 3 – maximum heating power → H1 LEDs are permanently lit
- Push 4 times → infrared program is activated → H1 LED double flashing interval


Infrared program:

In the infrared program, the infrared runs repeatedly through a wave from level 1 to level 3 during the previously set Timer time. In this way, an optimum infrared application in the different intensities is achieved.

When the time set in the Timer has elapsed, the heating elements are switched off in any mode – program or standby.

**Body zones**

---

## Pre-programmed intensity levels in the basic setting

Body zone	①	②	③	④	⑤
Button					
Area	Neck	Shoulder	Back/arms	Bottom/hands	Calves/legs
Intensity levels	<b>3</b>	<b>2</b>	<b>2</b>	<b>2</b>	<b>2</b>

## Programs P1 – P20

As a rule, the vibration units of the massage zones work with their pre-set intensity levels in all programs. The intensity of the vibration units, however, can of course be individually controlled or switched off in any program.

### Therapy programs according to Prof. Dr. Stutz

#### Basic programs

P01	Warm-up	For warming-up; loosening up muscular system.
P02	Deep Muscle vibration	Alleviating muscular pain; revitalising, enhancing general wellbeing.
P03	Lymphatic System Stimulation	Stimulating lymphatic flow; purification and deacidification as well as venous reflux.
P04	Pulse Vibration	Generates a stimulus that sets the body's fluids in motion and results in increased blood circulation in the affected tissue.
P05	Relaxation	Aimed at deep relaxation, stress reduction and revitalisation.
P06	Reflexology Stimulation	Stimulates blood and lymphatic circulation in case of circulatory disorders of calves and feet; pain reduction..

## Medical programs

P07	Tension Loosening	Treatment of back pain and tension of back muscles; treatment of myogelosis.
P08	Chronic Back Complaints	Aimed at chronic afflictions of lumbar spine as well as sciatic afflictions.
P09	Neck Relaxation	Apply for treatment of pain caused by neck or shoulder tension. Helps with shoulder pain..
P10	Joint Vitalisation	Aimed at joint pain, arthrosis, rheumatism and osteoporosis.
P11	Power Napping	Immediate relaxing effect.
P12	Activation	Setting the body's fluids in motion.
P13	Sports Recovery	After sports and physically exhausting activities.
P14	Specific Treatment of Visceral Fat	Additional use of tightened Andullation Belt – Transferring of Andullation Frequencies into abdominal area.
P15	Deep Relaxation	For stress reduction and relaxation; supporting meditation..
P16	Cardiovascular Stimulation	Stimulating vigour of the body
P17	Soothing Acute Pain	Apply for lumbago treatment or calf spasms, e.g.
P18	Sleep Disturbance Redressing	Aims for reduction of sleep disturbances; helps to reduce tension from everyday life.
P19	Body Tissue Relaxation	Removes blocks by reducing the tenison in the different layers of body tissue.
P20	Stress Tolerance	Reducing stress, Cool-Down.

## Safety instructions

---

The product complies with the basic requirements for medical devices according to the medical device directive 93/42/EEC.

**The manufacturer can, however, not be held liable if any damage or consequential damage is caused due to non-compliance with the following instructions:**

- Repair work may only be performed by the manufacturer. (see tech. data, page 5)
- The therapy system may only be used with the accessories which are included in the scope of delivery such as the infrared neck cushion, infrared back cushion, infrared calf cushion, andullation belt and power andullator. The indicated accessories may only be used in connection with the system. The accessories may not be used individually without the system or in a different way than intended.  
The use of different accessories and cables may have a negative influence on electromagnetic compatibility (EMC).
- All components and accessories of the therapy system are to be protected against moisture.
- Set the system up in a place where any electrostatic interactions with the device are avoided. Do not use it on plastic carpets or similar flooring!
- The area where the therapy system is used may not be left unattended since lines/cables pose a strangulation hazard for babies and children.
- Avoid tripping hazards caused by lines/cables and never support yourself on the manual control.
- If the treatment period is too long, contact injuries may result (redness, minor burns).
- The therapy system may only be transported in its original packaging.
- No changes may be made to the device which is not considered as intended use.
- The therapy system and all accessories must be protected against direct sunlight.
- Pets must be kept away from the device and its supply lines.
- The therapy system may not be used near heaters or fireplaces.
- If the product is not used as intended, is operated incorrectly or not properly repaired, the manufacturer will no longer be liable for any resulting damage.
- This product may only be used together with the power supply unit – see technical data, page 5.
- In case of errors during use, prior to cleaning and care and after use, always pull the mains plug of the power supply unit.
- Disconnect all poles of the device from the mains prior to cleaning, maintenance and disinfection.
- When you pull the mains plug out of the socket, never pull at the line and never use wet hands.
- Never pull the power supply line over sharp edges or jam it.
- Children and disabled people may not use the system without supervision.
- Do not use the product in a damp environment such as a bathroom or swimming pool.
- Do not stick sharp objects such as needles or similar in the device. Never transport the product in an unintended manner by lifting it at the connection lines.
- Never switch on the infrared heat under covers or cushions when the system is unattended.
- Depending on the frequency of use, the product is to be checked thoroughly for damage to the artificial leather cover and the cables from time to time.
- It is not allowed to disconnect the applied parts from the device during application. In addition, please ensure that the device is not accidentally disconnected during operation.

## User instructions

---

Regarding the storage and transport of the therapy system incl. accessories, the following ambient conditions are to be observed: (also applies to the unpackaged condition and interim storage)

Temperature range: -20 ... + 50 °C

Relative humidity: 15 ... 93 %, non-condensing

Atmospheric pressure: 700 ... 1060hPa (up to 3000 m)

Protect against moisture

### Caution:

If the location of the system is changed from a cold to a warm area, ensure that the device warms up slowly. To be on the safe side, only connect the system to the mains after 24 hours in the warm temperature area.

### Ambient conditions for the operation of the therapy system incl. accessories:

Temperature range: +5 ... + 40°C

Relative humidity: 15 ... 93 %, non-condensing

Atmospheric pressure 700 ... 1060hPa (up to 3000 m)

Protect against moisture

- For hygienic reasons, do not use the therapy system without your clothes on.
- You get a particularly intensive foot reflexology massage, if you additionally bend your legs and the soles of your feet are massaged by the vibration units of body zone 5.
- We recommend using the system for at least 15 minutes per application.
- There is no need to be worried if you are disturbed during your massage. When the fixedly programmed timer time has elapsed, your therapy system interrupts all set functions and switches to the safe standby mode.
- When the device is switched on or off between applications, it is immediately ready again for the intended use without any time delay.
- Please do not stand, sit or lie on the system when it is folded.
- Do not place any heavy or hot objects on the system.
- Extended direct exposure to sunlight or heating sources could lead to discolouration of the cover.
- If you remove the seal at the zipper, any warranty claims will become void.

### Contraindications

With the following contraindications, consult a physician before you use the therapy device: Persons with electronic implants, pregnant women and nursing mothers with thrombosis, inner lesions and bleeding, sensitivity disorders.

Please also observe the warnings regarding the use of infrared heat

## **Environmental protection**

---

Raw material recovery instead of waste disposal: device, accessories and packaging should be recycled in an environmentally compatible manner.

### Disposal information:

Proper disposal is environmentally friendly and avoids hazardous effects on people and the environment which may result from improper handling of waste products.

The packaging material is environmentally compatible and can be recycled. Please support us by disposing of the packaging in an environmentally compatible manner. The packaging may not be used as a children's toy. The folding cartons and films pose a suffocation hazard.

The device including accessories may not be disposed of in the household waste since it is made of high-quality materials which can be recycled and used again. According to the European directive 2012/19/EU (WEEE), electrical and electronic devices have to be collected separately and not together with the unsorted municipal waste so that they can be recycled.

Please obtain information on current ways of disposal from your city or municipal administration. The symbol with the crossed out wheelie bin indicates that separate collection is required.

Regarding the disposal of cleaning, care and disinfection agents, you must observe the information by the manufacturer of the product used.

Only applies to electronic components!

## **Power down**

---

- Press the OFF button: all functions are interrupted and the flashing LED next to the "OFF" button indicates readiness for operation.
- Now pull the mains plug from the socket.

## Cleaning, care and disinfection

---

Caution!

Pull the mains plug prior to cleaning, care or disinfection.

After cleaning or servicing, always clean or disinfect the device.

### Cleaning and care

Never immerse the device in water! The cover was made of a soft skai® material which is easy to clean. To ensure that you can enjoy your skai® material for a long time, please clean it regularly.

Soiling with oils, greases and seat has to be removed immediately. Use warm, mild suds and a microfibre cloth or a soft hand brush.

Do not use solvents, chlorides, polishing agents, chemical cleaning agents or wax polish.

Please note that staining by jeans or other textiles is not included in the manufacturer's warranty.

### Cleaning:

- Ensure that no residues of cleaning or care agents remain on the therapy system.
- Clean the therapy system after any soiling.

Caution: Do not use agents which contain oil, grease or alcohol!

This product is not resistant to dry cleaning.

Spots of ink have to be removed immediately.

### Disinfection:

- With clinical applications of the device or the accessories which require direct contact to the body, make sure that you disinfect the system prior to and after each use by wiping and observe the manufacturer's information on the disinfectant.
- Only the following agent is to be used for disinfection:  
Cleanisept Wipes 5% cloths by the company Schuhmacher GmbH.
- For disinfection/cleaning, please observe the data sheet/instructions by the manufacturer of the agent used.
- Use cloths by the company Schuhmacher GmbH to wipe all surfaces of the therapy system and the accessories. Only use cloths with an appropriate amount of moisture.

## Troubleshooting

The cause of an error does not necessarily have to be the product. You may save time and money, if you remedy simple errors yourself. The following information is supposed to help you with this:

SYMPTOM	CAUSE	MEASURES
<ul style="list-style-type: none"> <li>• OFF LED does not flash when device is ready for operation</li> </ul>	<ul style="list-style-type: none"> <li>▪ Mains plug of the power supply unit is not plugged in the socket</li> <li>▪ Massage is in sleep mode</li> <li>▪ Coupling connector of the manual control is not plugged in properly</li> <li>▪ Green LED at the power supply unit is not lit</li> <li>▪ Due to an error in the electronic control, the integrated fuse has tripped</li> <li>▪ Power supply unit is defective</li> <li>▪ Defective mains fuse or tripping of the automatic circuit breaker</li> <li>▪ Power outage</li> </ul>	<ul style="list-style-type: none"> <li>⇒ Plug in the mains plug</li> <li>⇒ Press any key</li> <li>⇒ Plug in the coupling connector</li> <li>⇒ Establish correct plug connection</li> <li>⇒ Repair by manufacturer</li> <li>⇒ See above</li> <li>⇒ Exchange mains fuse, or reset automatic circuit breaker</li> <li>⇒ -</li> </ul>
<ul style="list-style-type: none"> <li>• Massage function stops suddenly</li> </ul>	<ul style="list-style-type: none"> <li>▪ Timer time has elapsed</li> </ul>	<ul style="list-style-type: none"> <li>⇒ Restart</li> </ul>
<ul style="list-style-type: none"> <li>• Functions cannot be set or are not indicated by the corresponding LED</li> </ul>	<ul style="list-style-type: none"> <li>▪ Error of the electronic control</li> </ul>	<ul style="list-style-type: none"> <li>⇒ Repair by the manufacturer</li> </ul>
<ul style="list-style-type: none"> <li>• Vibration units are not functional</li> </ul>	<ul style="list-style-type: none"> <li>▪ Motor is jammed</li> </ul>	<ul style="list-style-type: none"> <li>⇒ See above</li> </ul>
<ul style="list-style-type: none"> <li>• Infrared heating cushions shut down prematurely</li> </ul>	<ul style="list-style-type: none"> <li>▪ Internal temperature regulation is activated, timer has elapsed</li> </ul>	<ul style="list-style-type: none"> <li>⇒ When the cushion has cooled down or the timer is re-activated, the function is restored</li> </ul>

## **Guidelines on electromagnetic compatibility – EMC guidelines**

---

This device was tested according to IEC 60601-1-2:2014 and meets all corresponding requirements.

Portable HF (high frequency) communication devices, including antennas, can have a negative influence on the device. Please operate the device (including accessories and cables) with a minimum distance of 30 cm from such HF devices.

The main performance characteristic of this electrical medical device is defined as application time with a maximum duration of 60 minutes. This main performance characteristic can cannot be influenced by the negative impact of electromagnetic radiation.

Regarding the basic and EMC safety, no special maintenance measures are required for this medical electrical device.

<p><b>The “VM 9100RM III therapy system” is intended for operation in the electromagnetic environment indicated below.</b></p> <p><b>The operator has to ensure that the “VM 9100RM III therapy system” is operated in this kind of environment.</b></p>		
<b>Measurements of the emitted electromagnetic interference</b>	<b>Conformity</b>	<b>Electromagnetic environment - Guidelines</b>
Radiofrequency emissions acc. to CISPR 11	Class B	<p><b>The “VM 9100RM III therapy system” is suitable for use at all facilities, including those in living areas and at facilities which are directly connected to a public supply network which also supplies buildings which are used for residential purposes.</b></p>
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Emission of voltage fluctuations and flickers acc. to IEC 61000-3-2	Compliant	

<p><b>The “VM 9100RM III therapy system” is intended for operation in the electromagnetic environment indicated below.</b></p> <p><b>The operator has to ensure that the “Therapy System VM 9100RM III” is operated in this kind of environment.</b></p>			
<b>Immunity testing</b>	<b>IEC 60601-1-2 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - Guidelines</b>
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	±8kV contact discharge ±15kV air discharge	±8kV contact discharge ±15kV air discharge	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity has to be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	±2kV for power cables ±1kV for input and output cables	±2kV for power cables ±1kV for input and output cables	The quality of the supply voltage should correspond to a typical business or hospital environment.
Surges acc. to IEC 61000-4-5	±1kV normal mode voltage ±2kV common mode voltage	±1kV normal mode voltage ±2kV common mode voltage	The quality of the supply voltage should correspond to a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations acc. to IEC 61000-4-11	<5% $U_T$ (>95% dip of $U_T$ ) for ½ period 40% $U_T$ (60% dip of $U_T$ ) for 5 periods 70% $U_T$ (30% dip of $U_T$ ) for 25 periods <5% $U_T$ (>95% dip of $U_T$ ) for 5 s	<5% $U_T$ (>95% dip of $U_T$ ) for ½ period 40% $U_T$ (60% dip of $U_T$ ) for 5 periods 70% $U_T$ (30% dip of $U_T$ ) for 25 periods <5% $U_T$ (>95% dip of $U_T$ ) for 5 s	If the user requires the continued functioning of the VM 9100RM III therapy system even if the energy supply is interrupted, it is recommended to supply the VM 9100RM III therapy system with power using an uninterruptible power supply or a battery.
Power frequency magnetic field (50-60Hz) acc. to IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields at mains frequency should correspond to the typical values as found in business and hospital environments.
<p>Comment: <math>U_T</math> = Alternating mains voltage prior to the application of the test levels.</p>			

**The “VM 9100RM III therapy system” is intended for operation in the electromagnetic environment indicated below.**  
**The operator has to ensure that the “VM 9100RM III therapy system” is operated in this kind of environment.**

Immunity testing	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - Guidelines
<p>Conducted radio-frequency acc. to IEC 61000-4-6</p> <p>Radiated radio-frequency acc. to IEC 61000-4-3</p>	<p>3 V<sub>eff</sub></p> <p>150 kHz to 80 MHz</p> <p>3 V/m</p> <p>80 MHz to 2.7 GHz</p>	<p>3 V<sub>eff</sub></p> <p>3 V<sub>eff</sub></p>	<p>Portable and mobile radio devices should not be used at a distance to the VM 9100RM III therapy system which is less than the recommended protective distance which is calculated based on the equation which is applicable for the transmission frequency.</p> <p><b>Recommended protective distance:</b></p> <p><math>d = 1.2 \text{ m } \sqrt{P}</math> for 80 MHz to 800 MHz</p> <p><math>d = 2.3 \text{ m } \sqrt{P}</math> for 800 MHz to 2.7 GHz</p> <p>P is the rated power of the transmitter in Watt (W) according to the information by the transmitter manufacturer and d is the recommended protective distance in metres (m). The field intensity of stationary radio transmitters should be less than the conformity level regarding all frequencies according to a test on site.</p> <p>In the surroundings of devices which bear the following symbol, interference is possible.</p> 

Comment 1: For 80 MHz and 800 MHz, the higher frequency range applies.  
 Comment 2: These guidelines may not be applicable in all cases. The spreading of electromagnetic quantities is influenced by the absorption and reflexions of buildings, objects and persons.

- a The field intensity of stationary transmitters such as base stations of radio telephones and mobile land radio equipment, amateur radios, AM and FM radio and television stations, cannot be exactly predetermined in theory. To determine the electromagnetic environment regarding the stationary transmitters, a study of the electromagnetic phenomena of the place of use should be considered. If the measured field intensity at the place of use of the VM 9100RM III therapy system exceeds the above-mentioned conformity levels, the VM 9100RM III therapy system should be observed to prove its functionality. If unusual performance characteristics are observed, additional measures may be required, e. g. different positioning or a different place of use of the VM 9100RM III therapy system.
- b Above the frequency range of 150 kHz to 80 MHz the field intensity should be below 3 V/m.

**Test specification for radiated radio-frequency IEC 61000-4-3:**

<b>Impulse frequency [MHz]</b>	<b>Band [MHz]</b>	<b>Service</b>	<b>Modulation</b>	<b>Max. power [W]</b>	<b>Distance [m]</b>	<b>Immunity test level [V/m]</b>
<b>385</b>	<b>380 – 390</b>	TETRA 400	Impulse modulation 18 Hz	<b>1.8</b>	<b>0.3</b>	<b>27</b>
<b>450</b>	<b>430 – 470</b>	GMRS 460, FRS 460	FM ± 5 kHz fluctuation 1 kHz sine	<b>2</b>	<b>0.3</b>	<b>28</b>
<b>710</b> <b>745</b> <b>780</b>	<b>704 – 787</b>	LTE Band 13, 17	Impulse modulation 217 Hz	<b>0.2</b>	<b>0.3</b>	<b>27</b>
<b>810</b> <b>870</b> <b>930</b>	<b>800 – 960</b>	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Impulse modulation  18 Hz	<b>2</b>	<b>0.3</b>	<b>27</b>
<b>1720</b> <b>1845</b> <b>1970</b>	<b>1700 – 1990</b>	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Impulse modulation  217 Hz	<b>2</b>	<b>0.3</b>	<b>27</b>
<b>2450</b>	<b>2400 – 2570</b>	Bluetooth, WLAN, 802,11 b/g/n, RFID 2450, LTE Band 7	Impulse modulation 217 Hz	<b>2</b>	<b>0.3</b>	<b>27</b>
<b>5240</b> <b>5500</b> <b>5785</b>	<b>5100 - 5800</b>	WLAN 802.11 a/n	Impulse modulation 217 Hz	<b>0.2</b>	<b>0.3</b>	<b>27</b>

<b>Recommended protective distances between portable and mobile radiofrequency telecommunication devices and the “VM 9100RM III therapy system”</b>			
The VM 9100RM III therapy system is intended for operation in an electromagnetic environment where the radiofrequency disturbances are controlled. The operators of the VM 9100RM III therapy system can help prevent electromagnetic interference by observing the minimum distances indicated below between portable and mobile radiofrequency telecommunication devices (transmitters) and the VM 9100RM III therapy system independent of the output power of the communication device.			
<b>Rated power of the transmitter</b>  (W)	<b>Protective distance, depending on the transmitter frequency (m)</b>		
	<b>150 kHz to 80 MHz</b>  $d= 1.2 \text{ m } \sqrt{P}$	<b>80 MHz to 800 MHz</b>  $d= 1.2 \text{ m } \sqrt{P}$	<b>800 MHz to 2.5 GHz</b>  $d= 2.3 \text{ m } \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters, the maximum rated power of which is not indicated in the table above, the recommended protective distance d in meters (m) can be determined by using the equation which belongs to the corresponding column where P is the maximum rated power of the transmitter in Watt (W) according to the information by the transmitter manufacturer.  Comment 1: For 80 MHz and 800 MHz, the higher frequency range applies. Comment 2: These guidelines may not be applicable in all cases. The spreading of electromagnetic quantities is influenced by the absorption and reflexions of buildings, objects and persons.			

**Service address**

Please send devices for service/repair to the following service address:

METEK GmbH  
Stammweg 8  
37327 Leinefelde-Worbis  
Phone: +49 (0)3605 5159-0  
Fax: +49 (0)36055159-59  
Email: [info@metek-massage.de](mailto:info@metek-massage.de)

**Manufacturer**

METEK GmbH  
Stammweg 8  
37327 Leinefelde-Worbis  
Germany

The medical device des is exclusively distributed under the brand name Andumedic® 3 by:

hhp  
home health products  
[www.hhp-international.com](http://www.hhp-international.com)

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